













<p><b>Step 1</b></p> <p>Complete date and time on tube and card Month in <u>letters</u> e.g. 01 JAN 2018</p> <div data-bbox="49 353 331 437"> <p>Trial Number: XXXX – XX Collection date: _/ _/ _ PROFILE – Unbuffered stool sample</p> </div>  	<p><b>Step 2</b></p>  <p>Please don't over-fill sample tube</p> 	<p><b>Step 3</b></p>  	<p><b>Step 4</b></p> 	<p><b>Step 5</b></p> <p><b>Important:</b> Check all contents are fully sealed and inside before closing. Once closed, cannot reopen.</p> 
<p><b>Step 6</b></p>  <p>Remove cardboard</p>	<p><b>Step 7</b></p> 	<p><b>Step 8</b></p> 	<p><b>Step 9</b></p> <p>Ready to be posted</p> 	<p><b>Step 10</b></p> <p>You can put First Class Safebox™ in any post-box.</p> 